

# A Robust, High Performance Revolutionary Demountable ICP Torch

- More robust and cost effective compared to other demountable torches.
- Demountable quartz or ceramic outer tube.
- Ceramic outer tube provides greater cost savings, long-life, and robust plasma conditions for increased sensitivity.
- Interchangeable injectors - choose from quartz, alumina, sapphire, and platinum with a range of diameters for aqueous, organic, high salt or HF-containing samples.
- Easier maintenance with the ability to remove the injector and outer tube.
- Direct replacement for standard torch.

Scan the QR Code to watch the  
D-Torch Installation Guide:

