

Eluo Nebulizer Cleaning Tool

Cleaning your nebulizer is as easy as 1, 2, 3!

1. Fill with solvent



2. Insert nebulizer



3. Clean nebulizer



Particle build-up in a nebulizer capillary and tip causes sample flow to be constricted, thereby reducing nebulizer efficiency and performance. Now, blocked nebulizers can be safely and easily restored to optimum performance with the Eluo.

The Eluo is designed to efficiently deliver a cleaning solution through the nebulizer capillary to dislodge particle build-up and thoroughly clean the nebulizer.

Scan the QR Code to watch the Nebulizer Maintenance video:



GLASS EXPANSION
Quality By Design