


Helix CT

ICP Spray Chamber

Consistent day-to-day analytical performance is just a click away with ConstantTorque™ technology.



The Helix CT spray chamber with ConstantTorque™ technology, provides a constant, reproducible, inert, gas-tight seal between the nebulizer and spray chamber. The Helix CT locking screw is equipped with a built-in torque control mechanism that allows for a consistent seal of the PTFE ferrule against the nebulizer, making it impossible to overtighten or under tighten while ensuring a gas-tight seal every time.

Scan the QR Code to watch the Helix CT Assembly video:



GLASS EXPANSION
Quality By Design