

Helix*

Распылительная камера без уплотнительных прокладок

Новинка!

- Полностью инертная поверхность камеры, устойчивая к воздействию органических растворителей
- Простота в установке и извлечении распылителя
- Камера имеет систему фиксации положения распылителя, гарантирующей его постоянство всякий раз при замене или настройке
- Сокращенное время промывки камеры благодаря отсутствию «мертвого» объема в зоне уплотнения распылителя



Даже наиболее устойчивые к химическому воздействию прокладки подвержены износу и разрывам. Они также имеют тенденцию «прилипать» к стеклянным распылителям, существенно усложняя процедуру их извлечения. Вот почему компания Glass Expansion разработала Helix – простую и удобную в обращении распылительную камеру без уплотнительных прокладок.

*Ожидает патента.



GLASS EXPANSION

Quality By Design

Head Office (Australia)

6 Central Boulevard
Port Melbourne VIC 3207
Australia
A.B.N. 35 006 560 032

Telephone: +61 3 9320 1111
Toll Free (Aust): 1800 777 638
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com
www.geicp.com

USA Office

4 Barlows Landing Road
Unit #2
Pocasset, MA 02559
USA

Telephone: 508 563 1800

Toll Free: 800 208 0097
Facsimile: 508 563 1802
Email: geusa@geicp.com

For your nearest Glass Expansion distributor, visit our website at www.geicp.com